Datasheet

Product Name: waveguard original cap, 128 channels with 5 HPI coil openings, 10/5, unshielded, Redel connector

Product Code: CA-069.s1

Product Version: 3

waveguard™ original 

Manufacturer:

eemagine Medical Imaging Solutions
Gubener Str. 47
10243 Berlin
Germany

Phone +49 (0)30 2904 8404

Fax +49 (0)30 2904 8405

E-Mail support@eemagine.com

Web www.eemagine.com


OQC passed

for part(s) with SERIAL number as indicated on product

Compatibility: MEGIN TRIUX MEG system

Available sizes: s1: S, M, L

Cap Layout:

Pinning:

Connector 1
yellow bend relief

Channel Number	Electrode	Connector Pin
1	Fp1	1
2	Fp2	2
3	F7	3
4	F3	4
5	Fz	5
6	F4	6
7	F8	7
8	T7	8
9	C3	9
10	Cz	10
11	C4	11
12	T8	12
13	P7	13
14	P3	14
15	Pz	15
16	P4	16
17	P8	17
18	O1	18
19	O2	19
20	Fpz	20
21	FCS	21
22	FC1	22
23	FC2	23
24	FC6	24
25	CP5	25
26	CP1	26
27	CP2	27
28	CP6	28
29	PO5	29
30	POz	30
31	PO6	31
32	Oz	32

Connector 2
green bend relief

Channel Number	Electrode	Connector Pin
33	AF7	1
34	AF3	2
35	AFz	3
36	AF4	4
37	AF8	5
38	F9	6
39	F1	7
40	F2	8
41	F10	9
42	FC3	10
43	FC4	11
44	C5	12
45	C1	13
46	C2	14
47	C6	15
48	CP3	16
49	CPz	17
50	CP4	18
51	P5	19
52	P1	20
53	P2	21
54	P6	22
55	PO3	23
56	PO4	24
57	FT7	25
58	FT8	26
59	TP9	27
60	TP7	28
61	TP8	29
62	TP10	30
63	PO7	31
64	PO8	32

Pinning:

Connector 3
blue bend relief

Channel Number	Electrode	Connector Pin
65	FT9	1
66	FT10	2
67	TPP9h	3
68	TPP10h	4
69	PO9	5
70	PO10	6
71	P9	7
72	P10	8
73	AFF1	9
74	AFF2	10
75	FFC5h	11
76	FFC3h	12
77	FFC4h	13
78	FFC6h	14
79	FCC5h	15
80	FCC3h	16
81	FCC4h	17
82	FCC6h	18
83	CCP5h	19
84	CCP3h	20
85	CCP4h	21
86	CCP6h	22
87	CPP5h	23
88	CPP3h	24
89	CPP4h	25
90	CPP6h	26
91	PPO1	27
92	PPO2	28
93	I1	29
94	I2	30
95	AFp3h	31
96	AFp4h	32

Connector 4
grey bend relief

Channel Number	Electrode	Connector Pin
97	M1	1
98	M2	2
99	AFF5h	3
100	AFF6h	4
101	FFT7h	5
102	FFC1h	6
103	FFC2h	7
104	FFT8h	8
105	FFT9h	9
106	FFT7h	10
107	FCC1h	11
108	FCC2h	12
109	FFT8h	13
110	FFT10h	14
111	TTP7h	15
112	CCP1h	16
113	CCP2h	17
114	FFT9h	18
115	FFT10h	19
116	CPP1h	20
117	CPP2h	21
118	TTP8h	22
119	PPO9h	23
120	PPO5h	24
121	PPO6h	25
122	PPO10h	26
123	POO9h	27
124	POO3h	28
125	POO4h	29
126	POO10h	30
127	OI1h	31
128	OI2h	32

Cable	
Cable length	1.8m
Cable	unshielded
Connector	
Connector type	Redel (32 pins)
Connector gender	Male

Warning:

Proper use of the CA-069 depends on careful reading of all instructions including the CA-XXX waveguard caps descriptions and labels that come with or on the devices. Inaccurate measurements may be caused by incorrect use of the device. Non-compliance with warnings and safety regulations may result in severe personal injury and total loss of equipment.

Disclaimer: We have attempted to write this document as accurately as possible. However, mistakes are bound to occur, and we reserve the right to make changes to the products, which may render parts of this document invalid. No part of this document may be copied or reproduced without the explicit permission of the authors.

English: Information for waveguard™ EEG cap and accessories

Manufacturer: eemagine Medical Imaging Solutions GmbH, Gubener Str. 47, 10243 Berlin, Germany

The **waveguard** EEG cap is a medical device of CF class I according to EU regulation 93/42/EEC, and also according to 2011/65/EU (RoHS). In the USA, **waveguard** EEG cap is a class II device according to the CFR, FDA 510(k) No. K110223. In Canada, certain cap types are registered as medical device class II according to the IMDR under MDL no. 88779.

Rx only

For product registrations in other countries send your inquiry to sales@eemagine.com.

Safety information

Intended use

- **waveguard** EEG caps are intended for recording of electroencephalographic (EEG) signals, including event-related potentials (ERP) and evoked potentials (EP) on subjects or patients with intact scalp in combination with medical-grade EEG recording equipment (CE class I/IIa certified inside the European Union).
- The cap must be dry and clean when used. The application of caps on subjects or patients should only be carried out by medically trained personnel (EEG technician). The technician must remain present during the entire recording session to prevent patient risks such as strangulation. Special infection control requirements may apply.

Contra-indications

- The caps should not be used on patients with skin injuries or problems such as blisters, burns, wounds from operation or other superficial wounds or transmissible diseases, such as – but not limited to – hepatitis, human immunodeficiency virus (HIV) or Creutzfeldt-Jacob disease.
- Do not use the cap in a Magnetic Resonance Imaging (MRI) scanner or in any other high-intensity electromagnetic field radiation environment unless the type of the cap expressly states the compatibility with such a device.
- In the event of defibrillation, disconnect the cap and/or electrodes from the EEG recording equipment first.

Side effects

- No side effects are known if properly applied. Please apply the correct size of the cap. Keep in mind that a too small cap may cause headache.
- Slight red marks on the exposed skin (e.g. at the forehead) due to contact pressure are a regular effect of wearing the cap and disappear within short time.

Important technical notes

- **Before the first use, make sure that the cap is properly disinfected and clean.**
- Do not use incompatible electrode materials like platinum or stainless steel in one recording.
- Do not use hard brushes or any other hard materials (in particular no metal tools such as screw drivers) to clean electrodes.
- Apply always the correct size of cap (related to head circumference).
- Make sure that the cap is completely dry before use.
- Make sure to use the appropriate amount of gel for each electrode; too much gel can lead to invalid signals.
- Do not clean the cap using an ultrasonic cleaning device.
- Do not keep the cap submerged in water for prolonged periods of time (maximum one to two hours).
- Do not expose the cap to hot environments over 50° Celsius (or short rinsing at 60° Celsius). The cap fabrics as well as the wiring do not withstand higher temperatures. Hot-air drying may significantly reduce the lifetime of the cap.
- Do not try to wring or squeeze the cap! We recommend to gently dab the cap fabric using a towel, or to hold the cap between two towels for drying.
- Do not use any autoclaving procedure on the cap.
- Make sure that the cap connectors always remain dry.
- Please dry the cap ideally in an upside-down position, with the cable being placed higher than the wet cap, in order to prevent adapter and tubing to be contaminated with water.
- Always disconnect caps holding the connector housing. Never drag at the cable! Note that dragged-out cables are not covered by the cap warranty.
- The sintered Ag/AgCl electrode elements do not need any rechlorification.

The User Guide contains additional information.

See www.ant-neuro.com/support or contact your local distributor to obtain a copy.

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